

Falcon® Cell Strainers and Containers



Falcon Cell Strainers

- ▶ A faster and easier alternative to gauze filtration in procedures involving dissociation of cells from either clumps or primary tissues
- ▶ Consistently obtain a more uniform single-cell suspension
- ▶ Three nylon mesh pore sizes for optimal performance in a variety of applications
- ▶ Sterilized and conveniently accessible in individual packaging
- ▶ Extended lip on strainer enables aseptic handling with forceps
- ▶ Strainers are made of a strong nylon mesh with 40, 70, or 100 µm pores that are evenly spaced for consistent results
- ▶ Molded color-coded polypropylene frame with tab enables easy handling

Features

- ▶ Sterilized by gamma irradiation
- ▶ Individually packaged

Falcon Cell Strainers Ordering Information

Cat. No.	Description	Qty/Pk	Qty/Cs
Falcon Cell Strainers			
352340	40 µm/Blue	1	50
352350	70 µm/White	1	50
352360	100 µm/Yellow	1	50

Tips

- ▶ A special sterile cell strainer assembly designed for cell sorting devices is listed on page 44. It contains a 35-micron nylon mesh strainer that is an integral part of the cap to a 5 mL polystyrene tube.



Falcon Containers

- ▶ Sterile, disposable polypropylene containers provide secure sample containment
- ▶ Available in 4½ oz. (110 mL) and 8 oz. (220 mL) sizes both with and without matching lids
- ▶ Feature molded-in graduations for easy measurements
- ▶ Graduated in ounces and milliliters to provide convenient, single-use, timesaving containers for collection, transportation, and storage of a wide variety of specimens
- ▶ Inert and chemically resistant to commonly used laboratory reagents at room temperature

Features

- ▶ Sterilized by gamma irradiation
- ▶ 4¼ oz. graduations from ½ oz. to 4½ oz. in ¼ oz. increments and from 20 mL to 110 mL in 10 mL increments
- ▶ 8 oz. size graduated from ¼ oz. to 8 oz. in ¼ oz. increments and from 20 mL to 220 mL in 10 mL increments

Falcon Containers Ordering Information

Cat. No.	Description	Qty/Pk	Qty/Cs
Falcon Containers			
354013	4½ oz. (110 mL) with Lid	1	100
354014	4½ oz. (110 mL) without Lid	20	500
354015	8 oz. (220 mL) with Lid	1	100
354020	8 oz. (220 mL) without Lid	20	500
354017	Lid for both sizes	20	500

Related Products

- ▶ Falcon Round-Bottom Tube with Cell Strainer Cap 56
- ▶ Falcon 50 mL Conical Tubes. 58