


SPECIAL OFFER

Kat.broj	Naziv	Jed.mj.	Nova cijena
638947	In-Fusion Snap Assembly MMix	10 Rxns	175 EUR-a
638948	In-Fusion Snap Assembly MMix	50 Rxns	700 EUR-a


- 1 Amplify**
each fragment/insert as a PCR product with 15-bp extensions (or 20-bp extensions for multiple fragment cloning), complementary to the linearized vector ends.



- 2 Combine**
in a single tube, the In-Fusion enzyme mix, the linearized vector, and PCR insert(s), and incubate for 15 min (liquid enzyme) or 30 min (EcoDry enzyme) at 50°C.



- 3 Transform**
competent cells with the In-Fusion cloning reaction. Plate on selective media and screen for positive clones.



MAKE CLONING EASY WITH IN-FUSION®



- cijene su bez pdv-a
- popusti ne vrijede za kupce koji imaju ugovorene uvjete
- cijene vrijede do kraja godine ili do isteka zaliha

Tel: +385 1 3369 610

+385 1 3369 611

E-mail: info@gorea-plus.hr